



LA TROBE
UNIVERSITY

MED3PRC

Practical Skills in Biomedical Science

Laboratory Workbook 2019

(Part 1: Pracs 1-3)

Name:

Student ID:

Total marks = / 93.5

Calculations

Dilutions

1. What is the dilution factor when 0.4 ml is added to 7.6 ml diluent? (1 mark)

2. How would you prepare 20 ml of a 1:400 dilution? (1 mark)

3. You have to make up sufficient antibody solution to add 100 μ l per well to a 96 well microtitre plate.

- a. How much antibody solution should you make (account for 5% overage)? (0.5 marks)

- b. How much stock antibody should you add to obtain a dilution of 1/20,000? (1 mark)

4. You have 0.6 mL of sample, and want to dilute it all to a fiftieth of its present concentration. How much diluent will you add, and what will the final volume be? (1 mark)

Solutions

5. You have 250 ml of 50X TAE stock solution. What volume of the 50X TAE stock is needed to prepare 35 ml of 10X TAE? (1 mark)




6. How would you prepare 800 ml of a 0.44 M NaCl solution (MW = 58.44 g/mole)? (1 mark)

7. You are asked to prepare HEPES buffered saline with the following constituents:

| | |
|-------------------|--------|
| HEPES | 10 mM |
| KCl | 5 mM |
| NaCl | 135 mM |
| CaCl ₂ | 1 mM |
| MgSO ₄ | 0.6 mM |
| Glucose | 6 mM |

You have previously made up **100 mM** stock solutions of HEPES, KCl, CaCl₂, and MgSO₄ and a **1 M** stock of NaCl. The MW of glucose is 180.16 g/mole

Calculate the volume of each solution, and the weight of glucose required to prepare 250 ml, 400 ml and 1 L of buffer. Fill in your answers in the table below (remember your volumes should be obtainable using standard laboratory equipment). (3.5 marks)

| | [final] (mM) | [Stock] (mM) | Volume (ml) | Volume (ml) | Volume (ml) |
|-------------------|-----------------|-----------------|---|---|---|
| HEPES | | | | | |
| KCl | | | | | |
| NaCl | | | | | |
| CaCl ₂ | | | | | |
| MgSO ₄ | | | | | |
| Glucose | | MW=180.16 g/mol |  g |  g |  g |
| H ₂ O | | | | | |
| Final volume | | | 250 | 400 | 1000 |

Protocols and Methods

- Using the recipe-like protocol provided for you in Pracs 4 and 5, construct a scientific method for publication by filling in the gaps below. (6 marks)

ELISA

_____ microliters of purified human PSA peptide (0-1000 ng/ml; Abcam AB41241) diluted in PBS (NaCl 137 mM; KCl 2.68 mM; Na₂HPO₄ 8 mM; KH₂PO₄ 1.5 mM; pH 7.3) was adsorbed onto microtitre plates for _____ at _____ °C. The plates were then washed _____ with _____ and blocked with _____ (w/v) low fat milk powder in PBS-T (PBS with _____ Tween 20) for _____ at _____ °C. Plates were washed _____ with _____ and incubated with _____ polyclonal anti-PSA peptide antibody Kallikrein loop (_____; Abcam® AB28563) for _____ at _____ °C. Wells were washed _____ with _____ and incubated with HRP-conjugated _____ anti-_____ IgG (_____; Sigma-Aldrich) for _____ at _____ °C. Wells were washed _____ with PBS (3x) and then developed using 1-step™ Ultra TMB-ELISA (Thermo Scientific). The colour reaction was stopped with _____ HCl and the final product read on a Spectrostar Nano plate reader at an optical density of _____.

Practical 1 – Expression of Bcl-2

AIM - PCR

Briefly state the aim of the PCR experiment (1 mark)

METHODS – PCR

Briefly outline the key steps in the method of the PCR experiment (see practicals 1 & 2) (3.5 marks)

EXPECTED RESULTS – PCR

Briefly indicate what you expect your results to be for each PCR reaction mix if the experiment was successful (2 marks)

| PCR reaction | Expected result |
|----------------------------|-----------------|
| No template control | |
| Positive control | |
| HL-60 ^{WT} DNA | |
| HL-60 ^{Bcl-2} DNA | |

LAB MANUAL CALCULATIONS

ASSESSING BCL-2 EXPRESSION AT THE GENE LEVEL

Part B - Quantify DNA

After extracting DNA from your cell pellet you are asked to determine the concentration of DNA in your lysate using the nanodrop. Fill out the table below to obtain 100 ng of DNA from your lysate. (1 mark)

| Cell culture (HL-60 ^{WT} or HL-60 ^{Bcl-2}) | DNA concentration (ng/μl) | Volume required for 100 ng of DNA (μl) |
|--|------------------------------|--|
| | | |

Part C - Prepare and run PCR reactions

Fill out the following table of reagents to make up the PCR reaction tubes. (2 marks)

(Reminder: you want to achieve a final concentration of 0.5 μM for each primer)

| Reagent | Volume (μl) | | |
|---|-------------|------------|------------|
| | Reaction 1 | Reaction 2 | Reaction 3 |
| Genomic DNA sample (~100 ng) or plasmid DNA control (~5ng) | | | |
| Sterile H ₂ O | | | |
| GoTaq™ Green Mix | 25 | 25 | 25 |
| Forward primer (5 μM stock) | | | |
| Reverse primer (5 μM stock) | | | |
| Total reaction volume: | 50 | 50 | 50 |

AIM – WESTERN BLOT

Briefly state the aim of the Western Blot experiment (1 mark)

METHODS – WESTERN BLOT

Briefly outline the key steps in the method of the Western Blot experiment (see practicals 1 & 2) (3.5 marks)

EXPECTED RESULTS – WESTERN BLOT

Briefly indicate what you expect your results to be for each sample tested if the Western Blot experiment was successful (3 marks)

| Primary antibody | Sample | | |
|------------------|----------------------------|-------------------------------|------------------|
| | HL-60 ^{WT} lysate | HL-60 ^{Bcl-2} lysate | Positive control |
| Anti-Bcl-2 | | | |
| Anti-FLAG | | | |

PRACTICAL 1 QUESTIONS

1. What is PCR? Briefly explain how it works considering the steps of denaturation, annealing and extension. (2 marks)

2. What is the purpose of the control reactions in the PCR – no template control and positive Bcl-2 control? (1 mark)

3. Approximately what size band would you expect to see if your PCR successfully amplified full length Bcl-2 gene? Why are we not amplifying the entire gene? (1 mark)

4. How can we tell the difference or detect transgenic Bcl-2 compared with a wildtype Bcl-2 gene by PCR? (1 mark)

Total marks for Practical 1 (Expression of Bcl-2):

/22 marks

Practical 2 – Expression of Bcl-2 (cont.)

LAB MANUAL CALCULATIONS

Part A – Prepare an agarose gel for PCR

For these calculations, refer to Prac 2 in your lab manual.

Indicate the amount of agarose you needed to weigh out to obtain a 0.8 % (w/v) solution when diluted with 0.5 X TBE. (1 mark)

Amount of agarose to obtain 0.8 % (w/v) solution: _____grams

Indicate the volume of SYBRSafe stain you need to add to your agarose/TBE solution to obtain 0.01 % (v/v). (1 mark)

Volume of SYBRSafe stain to obtain 0.01 % (v/v) solution: _____microliters

PRACTICAL 2 QUESTIONS

PCR

Refer to your DNA gel (or a gel deemed successful by the demonstrator) for the following questions.

1. Which samples(s) was there no PCR product? Why? (1 mark)

2. Why did we perform this PCR? (1 mark)

3. What is the difference between what we are detecting in the PCR and western blot today?
Why do we need to detect the presence of both? (1 mark)

4. If there was a PCR product for both HL-60^{WT} and HL-60^{Bcl-2} samples, where would the primers most likely bind to? (1 mark)

5. What is the overall conclusion from this PCR experiment? (1 mark)

Western Blot

6. What is the purpose of blocking? (1 mark)

7. Why is the secondary antibody anti-mouse? If a third antibody was needed to detect the secondary, what species would it be derived from? (1 mark)

8. How large is the Bcl-2 protein? What size band (in both HL-60^{WT} and HL-60^{Bcl-2} cells) are we looking for on the blot if probed with; anti-Bcl-2 or anti-FLAG antibodies? (2 marks)

| | HL60 WT | HL60 Bcl-2 |
|------------|---------|------------|
| Anti-Bcl-2 | | |
| Anti-FLAG | | |

Refer to your blot (or a blot deemed successful by the demonstrator) for the following questions. You should consider both anti-Bcl-2 and anti-FLAG probed blots.

9. Which sample(s) bound the anti-Bcl-2 antibodies? (1 mark)

10. Did the antibodies detect a polypeptide with the expected size for Bcl-2? How do you know this? (1 mark)

11. What is it about probing with either anti-Bcl-2 or anti-FLAG antibodies that helps distinguish between wildtype and transgenic Bcl-2? (1 mark)

12. Did we use an appropriate control for the anti-FLAG blot? If not, what should we use next time? (1 mark)

13. What is the overall conclusion from this Western Blot experiment? (1 mark)

Practical 3 – Caspase activity assay

AIM – CASPASE ACTIVITY ASSAY

Briefly state the aim of conducting the caspase activity assay in this experiment (1 mark)

METHODS – CASPASE ACTIVITY ASSAY

Briefly outline the key steps in the method of the caspase activity assay (3 marks)

EXPECTED RESULTS – CASPASE ACTIVITY ASSAY

Briefly indicate what you expect your results to be for each sample tested if the caspase activity assay experiment was successful (3 marks)

| Condition | Sample | | |
|----------------------------------|----------------------------|-------------------------------|------------------|
| | HL-60 ^{WT} lysate | HL-60 ^{Bcl-2} lysate | Negative control |
| Untreated | | | |
| Treated (with doxorubicin) | | | |

AIM – TRYPAN BLUE EXCLUSION

Briefly state the aim of conducting the caspase activity assay in this experiment (1 mark)

METHODS – TRYPAN BLUE EXCLUSION

Briefly outline the key steps in the method of the caspase activity assay (1.5 marks)

EXPECTED RESULTS – TRYPAN BLUE EXCLUSION


Briefly indicate what you expect your results to be (in terms of cell viability) for each sample tested if the trypan blue exclusion assay experiment was successful (2 marks)

| Condition | Sample | |
|----------------------------|----------------------------|-------------------------------|
| | HL-60 ^{WT} lysate | HL-60 ^{Bcl-2} lysate |
| Untreated | | |
| Treated (with doxorubicin) | | |

LAB MANUAL CALCULATIONS

Part A – preparing a mastermix

Complete the protocol table to prepare 500 µl of mastermix. (2 marks)

| Component | Stock concentration | Final conc in mix | Vol (µl) to add to final vol of 500 µl |
|----------------------------|---------------------|-------------------|--|
| HEPES buffer pH 7 | 1M | 0.1M | |
| PEG3350 | 50% | 10% | |
| CHAPS | 1% | 0.1% | |
| Dithiothreitol (DTT) | 1 M | 10 mM | |
| AcDEVD-AFC* | 2 mM | 50 µM | <i>*Do not add this but allow for its addition in your 500µL*</i> |
| Distilled H ₂ O | | |  |

DATA – CASPASE ACTIVITY ASSAY

3 marks

Refer to your data (or data deemed successful by the demonstrator) for the following questions.

Calculate the activity of caspase-3 as a change in fluorescence units/hr/1 x 10⁶ cells by following these steps.

- a. Calculate the change in fluorescence over the incubation period by subtracting the initial reading from the final reading (0.5 marks)
- b. Subtract the negative control value from sample values (0.5 marks)
- c. Calculate the average caspase-3 activity from your duplicate (1 mark)
- d. Repeat for the other cell type using another pair's initial and final readings. (1 mark)

MED3PRC Lab Workbook: Practical 3 – Caspase activity

| Cell type | Treatment | Initial Reading | Final Reading | Final - Initial | Final- initial-control | Duplicate average |
|------------------------|---------------|-----------------|---------------|-----------------|------------------------|-------------------|
| HL-60 ^{WT} | No treatment | | | | | |
| | | | | | | |
| | + doxorubicin | | | | | |
| | | | | | | |
| Negative control | | | | | | |
| HL-60 ^{Bcl-2} | No treatment | | | | | |
| | | | | | | |
| | + doxorubicin | | | | | |
| | | | | | | |
| Negative control | | | | | | |

DATA – TRYPAN BLUE EXCLUSION

4 marks

| Cell line = | | Cell number | | | | | %viability |
|-------------|----------------------|-------------|-------------|-----------|--------------|---------|------------|
| | Quadrant | Top Left | Bottom Left | Top Right | Bottom Right | Average | |
| Untreated | Trypan blue positive | | | | | | |
| | Trypan blue negative | | | | | | |
| Treated 24h | Trypan blue positive | | | | | | |
| | Trypan blue negative | | | | | | |
| Treated 48h | Trypan blue positive | | | | | | |
| | Trypan blue negative | | | | | | |


DATA COMPARISON

Obtain class data for trypan blue exclusion assay and determine for each condition (4 marks):

1. Overall average. In excel use the formula: =average(values)
2. Standard deviation. In excel use the formula: =stdev(values)

| | | HL-60 ^{WT} cells | | HL-60 ^{Bcl-2} cells | |
|-------------|----------------------|---------------------------|-------------------|------------------------------|-------------------|
| | | Class average | %viability +/- SD | Class average | %viability +/- SD |
| Untreated | Trypan blue positive | | | | |
| | Trypan blue negative | | | | |
| Treated 24h | Trypan blue positive | | | | |
| | Trypan blue negative | | | | |
| Treated 48h | Trypan blue positive | | | | |
| | Trypan blue negative | | | | |

Insert an appropriately labelled graph of YOUR caspase activity assay results (3 marks)



*To correctly insert your graph, copy your completed graph from Excel into Powerpoint then right click on graph and "Save Picture As" as JPEG file, then upload the JPEG file. It is your responsibility that the graph is correctly uploaded

Insert an appropriately labelled graph of the CLASS trypan blue exclusion assay results (with SD)
(3 marks)

PRACTICAL 3 QUESTIONS

1. Why does measurement of caspase-3 activity give you an indication of apoptotic cell death?
(1 mark)

2. Why were cells treated with doxorubicin? (1 mark)

3. Compare fluorescence of lysates from treated and untreated cells for each cell type. What does this assay tell you about doxorubicin treatment? (2 mark)

4. How does cell viability change with a longer incubation in doxorubicin (24h vs 48h)? What effect does Bcl-2 over expression have on this? (1 mark)

5. Compare your observations for the caspase activity assay and the trypan blue exclusion assay. Do both assays support the same trend in regards to Bcl-2 over expression? (1 mark)

6. Can we make conclusions on the effect of Bcl-2 over expression on apoptosis by only considering the trypan blue exclusion data? Why? (2 marks)

7. What is the overall conclusion from this caspase activity assay? (1 mark)

8. What is the overall conclusion from this trypan blue exclusion assay? (1 mark)

Total marks for Practical 3 (Caspase activity assay):

/39.5 marks

